

## APPLICATION NOTE

# Viral Clearance under Mild pH using Protein A Resins

Alvaro Cruz-Izquierdo (1), Mark Hicks, Greg Saunders

## Introduction

Host Protein A affinity resins are essential for purifying monoclonal antibodies (mAbs) and similar therapeutic molecules. During the purification process, the target molecule binds to the affinity column, and a lower pH is typically used for elution. However, using typical elution levels (pH 3–pH 3.5) can result in sensitivity issues, often manifested as increased aggregation or other unwanted modifications – significantly reducing yield.

## Study Overview

In this study, the performance of two Purolite™ chromatography resins, Praesto Jetted A50 and Praesto Jetted A50 HipH, was compared to that of a commercially available agarose-based protein A resin (Competitor A).

This study evaluated the viral clearance capabilities of Praesto Jetted A50 under standard protein A elution conditions, as well as Praesto Jetted A50 HipH at both standard and milder elution levels. The viruses and experimental conditions utilized in this study were chosen based on industry guidelines (ICH Q5A), with the goal of determining log reduction values (LRV) for viral clearance for two model viruses: one enveloped (MLV) and one non-enveloped (MMV).

A total of sixteen runs were performed, two per virus across each elution condition studied. Duplicate runs for each virus and condition were conducted to verify the consistency of the reduction factors. Consistency was shown by the log reduction factors of each duplicate run being within one log of each other.

**TABLE 1 Summary of Viral Clearance Runs**

Summary of runs performed and elution conditions.

| Resin                   | Elution | MLV      | MMV      |
|-------------------------|---------|----------|----------|
| Praesto Jetted A50      | pH 3.5  | 2        | 2        |
| Praesto Jetted A50 HipH | pH 3.5  | 2        | 2        |
| Praesto Jetted A50 HipH | pH 4.5  | 2        | 2        |
| Competitor A            | pH 3.5  | 2        | 2        |
| <b>TOTAL RUNS/VIRUS</b> |         | <b>8</b> | <b>8</b> |

**TABLE 2** Chromatography Column Parameters

| Parameter                | Small-Scale Model                                  |
|--------------------------|--|
| Column Diameter (cm)     | 1.0  |
| Column Bed Height (cm)   | 20   |
| Column Volume            | 15.7 mL  |
| Linear Flow Rate (cm/h)  | 300 for process<br>200 for loading                 |
| Volumetric Flow Rate     | 3.93 mL/min for process<br>2.62 mL/min for loading |
| Load Concentration (g/L) | 2.8  |
| Load Ratio               | 45 g of product per L of resin                     |

**TABLE 3** Small-Scale Method for Viral Clearance Study

| Step                 | Solution  | Target Volume (CV) | Flow Rate (cm/h) | Estimated Time (mins) |
|----------------------|---|--------------------|------------------|-----------------------|
| Rinse                | Purified Water  | 2                  | 160              | 15                    |
| Pre-Use Sanitization | 0.5 M Sodium Hydroxide  | 3                  | 240              | 15                    |
| Rinse                | Purified Water  | 2                  | 300              | 8                     |
| Neutralization       | 20 mM Sodium phosphate, 150 mM sodium chloride, pH 7.4            | 2                  | 300              | 8                     |
| Equilibration        | 20 mM Sodium phosphate, 150 mM sodium chloride, pH 7.4            | 4                  | 300              | 15                    |
| Load                 | Pro A Load<br>Target load ratio: 45 g/L                           | Variable           | 200              | Variable              |
| Post Load Wash 1     | 20 mM Sodium phosphate, 500 mM sodium chloride, pH 7.4            | 1.5                | 200              | 9                     |
|                      |   | 3.5                | 300              | 14                    |
| Post Load Wash 2     | 50 mM Tris, 1 M sodium chloride pH 8.0                            | 5                  | 300              | 20                    |
| Post Load Wash 3     | 50 mM Sodium phosphate, pH 7.0                                    | 3                  | 300              | 12                    |
| Elution              | 50 mM Sodium acetate, pH 3.5<br>Or<br>50 mM Sodium acetate pH 4.5 | 6                  | 300              | 24                    |
| Strip                | 100 mM Acetic Acid  | 2                  | 300              | 8                     |
| Rinse                | Purified Water  | 2                  | 300              | 12                    |

**TABLE 4 Total Virus Loaded for Small-Scale Model Runs**


| Resin                                    | MMV (by TCID50)   | MLV (by QPCR)      |
|--|-------------------|--------------------|
| Praesto Jetted A50                       | 8.34 ± 0.36 log10 | 11.35 ± 0.04 log10 |
| Praesto Jetted A50 HipH (Elution pH 3.5) | 8.43 ± 0.33 log10 | 11.40 ± 0.04 log10 |
| Praesto Jetted A50 HipH (Elution pH 3.5) | 8.25 ± 0.38 log10 | 11.39 ± 0.04 log10 |
| Competitor A                             | 8.52 ± 0.28 log10 | 11.23 ± 0.04 log10 |

## Viral Clearance Results

Overlays of the UV A280nm profile of the chromatographs for the MLV-spiked runs are shown in Figures 1–4 and MMV-spiked runs are shown in Figures 5–8. Each figure includes both runs performed at that elution condition (first in blue, second in red).


**Figure 1**

### MLV-Spiked Run for Praesto Jetted A50 at pH 3.5




**Figure 2**

### MLV-Spiked Run for Praesto Jetted A50 HipH at pH 3.5




**Figure 3**

### MLV-Spiked Run for Praesto Jetted A50 HipH at pH 4.5




**Figure 4**

### MLV-Spiked Run for Competitor A at pH 3.5




**Figure 5**

**MMV-Spiked Run for Praesto Jetted A50 at pH 3.5**




**Figure 8**

**MMV-Spiked Run for Competitor A at 3.5**



**Figure 6**

**MMV-Spiked Run for Praesto Jetted A50 HipH at pH 3.5**




**Summary**

Viral clearance results indicate significant clearance across each protein A chromatography resin studied – with all resins showing a log reduction factor greater or equal to 2 log<sub>10</sub>. Duplicate runs proved consistent as all duplicate runs remained within 1 log of each other. In this study, Praesto Jetted A50 HipH exhibited better removal of both viruses compared to Praesto Jetted A50 and the competitor resin. For MLV, the virus log reduction of Praesto Jetted A50 HipH with high elution pH 4.5 showed better removal (log reduction of 4.85 ± 0.06 log<sub>10</sub>) compared to any run using the typical elution of pH 3.5. For MMV, the virus log reduction of Praesto Jetted A50 HipH with high elution pH 4.5 showed better removal (log reduction of 3.85 ± 0.52 log<sub>10</sub>) compared to any run using the typical elution of pH 3.5.

**Figure 7**


**MMV-Spiked Run for Praesto Jetted A50 HipH at pH 4.5**



Ecolab is a global developer, manufacturer, and supplier of Purolite™ Resins including ion exchange, catalyst adsorbent and advanced polymers that make the world cleaner and healthier.



[www.puroliteresins.com](http://www.puroliteresins.com)



## We're ready to solve your process challenges.

For further information on products and services, visit [www.puroliteresins.com](http://www.puroliteresins.com) or contact us at the addresses below.

### Americas

[americas@ecolab.com](mailto:americas@ecolab.com)

### Asia Pacific

[asiapacific@ecolab.com](mailto:asiapacific@ecolab.com)

### EMEA

[emea@ecolab.com](mailto:emea@ecolab.com)

The statements, technical information and recommendations contained herein are believed to be accurate as of the date hereof. Since the conditions and methods of use of the product and of the information referred to herein are beyond our control, Purolite expressly disclaims any and all liability as to any results obtained or arising from any use of the product or reliance on such information; NO WARRANTY OF FITNESS FOR ANY PARTICULAR PURPOSE, WARRANTY OF MERCHANTABILITY OR ANY OTHER WARRANTY, EXPRESSED OR IMPLIED, IS MADE CONCERNING THE GOODS DESCRIBED OR THE INFORMATION PROVIDED HEREIN. The information provided herein relates only to the specific product designated and may not be applicable when such product is used in combination with other materials or in any process. Nothing contained herein constitutes a license to practice under any patent and it should not be construed as an inducement to infringe any patent and the user is advised to take appropriate steps to be sure that any proposed use of the product will not result in patent infringement.



©2024 Purolite  
All rights reserved.  
P-000152-100PP-92024-ENG-R1-BP